

Technology	Rapid PCR combined with proprietary electrochemical detection																			
Regulatory designation	FDA 510(k) cleared Moderately complex with point of care designation CLIA-Waiver anticipated soon																			
Targets	Chlamydia (CT) genomic DNA, single target and gonorrhea (NG) genomic DNA, dual target <ul style="list-style-type: none"> CT: evaluated against a panel of 15 serovars, including the new or Swedish variant (nvCT) NG: evaluated against a panel of 32 different strains including strains from WHO, geographically diverse clinical isolates and strains with known mutations 																			
Specimen types	<p>Female vaginal swab</p> <ul style="list-style-type: none"> Clinician-collected Self-collected by a patient in a clinical setting <p>Male first-catch urine</p> <p>All specimen types are cleared for asymptomatic and symptomatic patients</p>																			
Hands-on time	<1 minute																			
Precision pipetting	Not required																			
Turnaround time	~30 minutes																			
Results	Easy to understand qualitative results with no interpretation required (Detected/Not Detected output for CT/NG)																			
Clinical performance	<table border="1"> <thead> <tr> <th rowspan="2">Target</th> <th colspan="2">Female</th> <th colspan="2">Male</th> </tr> <tr> <th>Sensitivity</th> <th>Specificity</th> <th>Sensitivity</th> <th>Specificity</th> </tr> </thead> <tbody> <tr> <td>Chlamydia</td> <td>96.1%</td> <td>99.1%</td> <td>92.5%</td> <td>99.3%</td> </tr> <tr> <td>Gonorrhea</td> <td>100.0%</td> <td>99.9%</td> <td>97.3%</td> <td>100.0%</td> </tr> </tbody> </table>	Target	Female		Male		Sensitivity	Specificity	Sensitivity	Specificity	Chlamydia	96.1%	99.1%	92.5%	99.3%	Gonorrhea	100.0%	99.9%	97.3%	100.0%
Target	Female		Male																	
	Sensitivity	Specificity	Sensitivity	Specificity																
Chlamydia	96.1%	99.1%	92.5%	99.3%																
Gonorrhea	100.0%	99.9%	97.3%	100.0%																
Limit of Detection (LoD)	<table border="1"> <thead> <tr> <th rowspan="2">Organism</th> <th colspan="2">Genome equivalents/mL</th> </tr> <tr> <th>Vaginal swab specimens</th> <th>Male urine specimens</th> </tr> </thead> <tbody> <tr> <td>CT serovar E (ATCC-VR-348B)</td> <td>407.4</td> <td>484.3</td> </tr> <tr> <td>CT serovar F (ATCC-VR-346)</td> <td>755.5</td> <td>769.3</td> </tr> <tr> <td>NG strain ATCC 49226</td> <td>245.6</td> <td>125.6</td> </tr> <tr> <td>NG strain ATCC 700825</td> <td>206.1</td> <td>212.3</td> </tr> </tbody> </table>	Organism	Genome equivalents/mL		Vaginal swab specimens	Male urine specimens	CT serovar E (ATCC-VR-348B)	407.4	484.3	CT serovar F (ATCC-VR-346)	755.5	769.3	NG strain ATCC 49226	245.6	125.6	NG strain ATCC 700825	206.1	212.3		
Organism	Genome equivalents/mL																			
	Vaginal swab specimens	Male urine specimens																		
CT serovar E (ATCC-VR-348B)	407.4	484.3																		
CT serovar F (ATCC-VR-346)	755.5	769.3																		
NG strain ATCC 49226	245.6	125.6																		
NG strain ATCC 700825	206.1	212.3																		
External controls	Commercially available through ZeptoMetrix (NATtrol™ Cat Nos. NATCT(434)-6MC, NATNG-6MC)																			
Internal controls	Included																			
Preventative maintenance	Not required																			
Calibration	Not required																			
Platform specifications	0.9 ft x 0.9 ft x 1.2 ft (H x W x D) 22.5 lbs																			
Cartridge storage	Refrigerated cartridge storage; ambient cartridge storage coming soon																			
Instrument storage	Ambient temperature																			
CPT Codes	Chlamydia: 87491 - Gonorrhea: 87591																			
Tests in development, including:	<i>Trichomonas vaginalis</i> , <i>Mycoplasma genitalium</i> , and SARS-CoV-2																			

The binx health *io* CT/NG Assay, when tested using the binx health *io* Instrument, is a fully automated, rapid, qualitative test intended for use in point-of-care or clinical laboratory settings for the detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* DNA by polymerase chain reaction. The binx health *io* CT/NG Assay is intended for use with female vaginal swab specimens, collected either by a clinician or self-collected by a patient in a clinical setting, or male urine specimens, as an aid in the diagnosis of symptomatic or asymptomatic *Chlamydia trachomatis* and/or *Neisseria gonorrhoeae* infection. For a symptomatic male patient with a chlamydia negative test result, further testing with a laboratory-based molecular test is recommended.

510(k) clearance does not constitute approval by FDA of a device.