

Technology Ultra-Rapid PCR combined with proprietary electrochemical detection

Regulatory designation FDA 510(k) cleared
 CLIA-Waived

Targets Chlamydia (CT) genomic DNA, single target and gonorrhea (NG) genomic DNA, dual target

- CT: evaluated against a panel of 15 serovars, including the new or Swedish variant (nvCT)
- NG: evaluated against a panel of 32 different strains including strains from WHO, geographically diverse clinical isolates and strains with known mutations

Sample types

Female vaginal swab

- Clinician-collected
- Self-collected by a patient in a clinical setting

Male first-catch urine

Both sample types are cleared for asymptomatic and symptomatic patients

Hands-on time <1 minute

Precision pipetting Not required

Turnaround time ~30 minutes

Results Easy to understand qualitative results with no interpretation required
 (Detected/Not Detected output for CT/NG)

Clinical performance	Female		Male	
	Sensitivity	Specificity	Sensitivity	Specificity
Chlamydia	96.1%	99.1%	92.5%	99.3%
Gonorrhea	100.0%	99.9%	97.3%	100.0%

Limit of Detection (LoD)	Genome equivalents/mL		
	Organism	Vaginal swab specimens	Male urine specimens
CT serovar E (ATCC-VR-348B)		407.4	484.3
CT serovar F (ATCC-VR-346)		755.5	769.3
NG strain ATCC 49226		245.6	125.6
NG strain ATCC 700825		206.1	212.3

External controls Commercially available through ZeptoMetrix (NATtrol™ Cat. Nos. NATCT(434)-6MC, NATNG-6MC)

Internal controls Included

Preventative maintenance Not required

Calibration Not required

Platform specifications 0.9 ft x 0.9 ft x 1.2 ft (H x W x D)
 22.5 lbs

Cartridge storage Refrigerated cartridge storage; room temperature cartridge storage coming soon

Operation Temperatures Optimal operation temperatures are 10°C to 35°C

CPT Codes Chlamydia: 87491 - Gonorrhea: 87591

Tests in development, including: *Trichomonas vaginalis*, *Mycoplasma genitalium*, and SARS-CoV-2

The binx health io CT/NG Assay, when tested using the binx health io Instrument, is a fully automated, rapid, qualitative test intended for use in point-of-care or clinical laboratory settings for the detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* DNA by polymerase chain reaction. The binx health io CT/NG Assay is intended for use with female vaginal swab specimens, collected either by a clinician or self-collected by a patient in a clinical setting, or male urine specimens, as an aid in the diagnosis of symptomatic or asymptomatic *Chlamydia trachomatis* and/or *Neisseria gonorrhoeae* infection. For a symptomatic male patient with a chlamydia negative test result, further testing with a laboratory-based molecular test is recommended.

510(k) clearance does not constitute approval by FDA of a device.